
PORCINE CARDIOMYOCYTE ISOLATION SOLUTIONS COMPOSITION

Adam B. Veteto

IonOptix Protocols

09/30/2019



LIST OF CHEMICALS AND REAGENTS

Chemical	Company	Catalog #
Minimum Essential Medium Eagle (MEM)	Sigma	M0518-10X1L
NaHCO ₃	Sigma	S8875-500G
Pyruvic Acid	Sigma	P-2256
Na-HEPES	Sigma	H7006-500G
HEPES	Sigma	H3375-1KG
Heparin		
Penicillin-Streptomycin (10,000U pen + 10,000mg strep per mL)	Gibco/Invitrogen	15140-122
Albumin from bovine serum (BSA)	Sigma	A6003-100G
Taurine	Sigma	T0625-100G
Liberase TH Research Grade	Roche	54 011 510 001



SOLUTIONS TO PREPARE PRIOR TO DAY OF EXPERIMENTS

Liberase TH Research Grade aliquots

- Dissolve 50 mg of Liberase TH in 10 mL dH₂O
- Sit on ice to dissolve for 30 minutes
- Separate into 1.125 mL aliquots (2.25 mL of 5 mg/mL = 11.25 mg)
2.25 mL will be added to 500 mL enzyme perfusion solution
- Freeze at -80 until day of experiment

Zero Ca²⁺ PSS Base Solution

<u>Chemical</u>	<u>FW or fluid conc</u>	<u>Amount /2L</u>	<u>final [mM]</u>
NaCl	58.44	16.71 g	143
KCl	74.55	0.7455 g	5
d-Glucose	180	3.6 g	10
HEPES	238.3	4.766 g	10
MgCl ₂	1 M sln	2 mL	1

fill to ~2L
 stir, pH with NaOH to 7.35, fill to 2L
 Store at 4 degrees, make fresh weekly

MEM Base Solutions

<u>Chemical</u>	<u>FW or fluid conc</u>	<u>Amount/L</u>	<u>final [mM]</u>
MEM	-	11.19 g	1X
PenStrep sln	10,000 U/mL	5 mL	50,000 U/L
NaHCO ₃	84.01	0.84 g	10
Na-Pyruvate	110.04	0.22 g	2
Na-HEPES	260.29	2.6 g	10
HEPES	238.3	2.38 g	10

fill to ~1L
 stir, pH to 7.35, fill to 1L
 Filter sterilize with 0.22 µm SteriCups
 Store at 4 degrees, keep sterile, make fresh weekly



DAY-OF SOLUTIONS

PORCINE CARDIOMYOCYTE ISOLATION: LEFT-VENTRICULAR WEDGE PREPARATION

1

Solution 1: Initial Wash

To **500ml of 0 Ca²⁺ Base Solution**, add:

1 ml of heparin (1000 U/ml)

Final concentration of solution:

NaCl	143 mM
KCl	5 mM
d-Glucose	5 mM
Hepes	5 mM
MgCl ₂	1 mM
Heparin	2 U/mL

2

Solution 2: Enzyme Perfusion Solution

To **500 mL MEM Base solution** add:

500 mg Taurine (MW 125.1)

100 µl of 0.1 M CaCl₂ solution

1 ml of heparin (10,000U/10ml)

Heat to 37 degrees, just prior to use add 2.25 mL of Blendzyme TH solution.

Final concentration of solution:

MEM	1X
NaHCO ₃	7 mM
Na-Pyruvate	2 mM
Na-HEPES	7 mM
HEPES	7 mM
CaCl ₂	20 µM
Taurine	10 mM
PenStrep	50,000 U/L
Blendzyme TH	0.045 mg/mL

Sol 1:15 min
Sol 2: 28 min

3

Solution 3: BSA Solution

To **200 mL MEM Base solution** add:

2 g BSA

200 µl of 0.1 M CaCl₂ solution.

Final concentration of solution:

MEM	1X
NaHCO ₃	10 mM
Na-Pyruvate	2 mM
Na-HEPES	10 mM
HEPES	10 mM
CaCl ₂	50 µM
PenStrep	50,000 U/L
BSA	10 mg/mL

4

Solution 4: Final Wash Solution

To **100 mL MEM Base solution** add:


125 mg Taurine (MW 125.1)

48 µl of 0.1 M CaCl₂ solution.

Final concentration of solution:

MEM	1X
NaHCO ₃	10 mM
Na-Pyruvate	2 mM
Na-HEPES	10 mM
HEPES	10 mM
CaCl ₂	50 µM
Taurine	10 mM
PenStrep	50,000 U/L

LAD Cannulation



Coronary ligations

ONOPTIX

Myocyte Isolation


VENTRICULAR MYOCYTES (MID-MYOCARDIUM):

- At ~30 minutes post digestion, cut through center of LV wedge. If color (orange) and texture (soft but not slimy) of a region is appropriate, cut out 1 inch square section and place in Solution 3. Remove epi and endocardium, and place mid-myocardial region into fresh Solution 3 in a Medium weight boat. Cut into smaller segments with fine scissors. (Note: Sharpness of scissors matters, otherwise you crush cells and live/dead ratio will drop. Also, don't dice tissue, just cut into smaller segments so cells can slough off). Gently shake tissue with forceps for 3 minutes and cells will come off.

- Filter cells through 200 μ m Nylon Filter Mesh into 50 mL Falcon Tube
- Place Falcon Tube (lid on but open to atmosphere) at ~30 degree angle for 5 minutes, and cells will settle on the bottom/side of tube.
- Turn tube upright and gently twist tube by hand to release cells from the side, and they will gradually form a soft pellet at the bottom of the tube.
- Immediately Resuspend pellet in 35 mL *Solution 3*, and check density and quality of cells.
- Allow cells to settle (30 degrees) and pellet (twist tube) again.
- Resuspend this pellet in 10-40 mL of *Solution 4*, depending on how dense you want the cells.



Filter through 0.2 mm Nylon mesh into 50 mL Falcon Tube



Example ventricular myocyte pellet